# rinowash Ego

INSTRUCTIONS FOR USE • ISTRUZIONI PER L'USO  
INSTRUCTIONS POUR L'EMPLOI • ANWEISUNGEN FÜR DEN GEBRAUCH  
INSTRUCCIONES DE USO • INSTRUÇÕES DE USO • INSTRUKCJA OBSŁUGI

 **Air Liquide**  
HEALTHCARE

ENGLISH • ITALIANO • FRANÇAIS • DEUTSCH • ESPAÑOL • PORTUGUES • POLSKI • РУССКИЙ

Fig. A







Fig. B


Fig. C

Fig. D


Fig. E


Fig. F



Fig. G

# rinowashEgo

## INTENDED USE

The device is a micronized nasal irrigation device for aerosol therapy, which allows treating disorders of the upper respiratory tract and physiological washing (with saline, hypertonic solution or thermal water) and/or therapeutic washing (with medicines) of nasal cavities.

The device produces a micronized solution jet which, when directed into the nasal cavities, promotes hydration, fluidification and the removal of mucus or phlegm.

 *The administration of medicines should be subject to medical prescription which defines the type of medicine, the doses to be administered and the treatment duration.*

Rinowash Ego has been designed and manufactured according to the provisions laid down in the Directive 93/42 EEC on medical devices (with subsequent updates) and it bears the CE marking.

## DESCRIPTION OF SYMBOLS USED

	Class II
	Type BF applied part
<b>IP22</b>	Device protected against contact with fingers and drops of water with a maximum inclination of 15°.
	Correct disposal of the product: The product is compliant with the regulation on electrical and electronic equipment and should not be disposed together with domestic waste.
<b>REF</b>	Product code
<b>SN</b>	Serial Number
<b>CE</b> 0051	CE Marking
	Temperature limits

	Relative humidity limits
	Atmospheric pressure limits
	Read the instructions for use
	Manufacturer
	Keep dry
	Medical Device

## WARNINGS

- Carefully read the instructions for use and keep them in a safe place.
- Do not pour water directly onto the device and do not immerse it in water
- Do not use the device while taking a bath or shower.
- The device should be used only as described in this manual. Any other use is not permitted as is considered improper and dangerous.
- Some components, due to their small size, if swallowed, may pose a risk of suffocation. Always use the device under adult supervision.
- The power supply cable, due to its length, may cause strangulation. Always use the device under adult supervision.
- It is recommended to check that in the table (Chapter “Disposal of the device and its components”), there are no materials for which some form of allergic reaction has occurred in the past.
- The use of the device by a child should always occur under adult supervision who is aware of these instructions.
- The administration of medicines should be subject to medical prescription which defines the type of medicine, the doses to be administered and the treatment duration.
- Do not use essential oil solutions (for example, menthol, eucalyptol, etc.) as they are not compatible with the materials of the device.
- Do not use the device while taking a bath or shower.
- The manufacturer has defined a service life of at least 400 hours from first use.
- In case of failure and/or malfunctioning of the device, please refer to the chapter “Troubleshooting”. Do not tamper with or open the base unit of the device. For repair operations, contact only a technical service centre

authorised by the manufacturer. Failure to observe the above provisions may compromise the safety of the device.

- In case of doubts on the operation, if you need assistance, or you want to report unforeseen events, contact the manufacturer.
- The device contains the battery and must be kept away from heat sources, please read the conditions of use and storage.
- Do not touch the electrical outlet and electrical components of the device with wet hands.
- Device protected against access to hazardous parts with a finger and against vertically falling drops of water (IP22).
- In case of an accidental fall into water, the device can be extracted only after interrupting the power supply. After such an event, the device may not be used and requires a complete overhaul.
- The compressed air generated and available at the air outlet connection can be a potential hazard.
- The device is not suitable for use in anaesthesia and pulmonary ventilation systems.
- The device is not suitable for use in the presence of anaesthetic mixtures with air, oxygen or nitrous oxide.
- The performance of the device is guaranteed if the accessories used are the original ones, always use original accessories. For safety reasons, use original USB cable and power supply adapter.
- The correct operation of the device may be compromised by electromagnetic interferences exceeding the limits expressed by the standards in force. If the device interferes with other electrical devices, remove the device.
- Do not pull the power supply cable to unplug it from the wall outlet and do not place the device in a location where it will be difficult to unplug.
- The device only operates on battery power, it is not possible to use the device when it is connected to the power supply network for charging.
- The battery is not replaceable.
- Once charging is complete, disconnect the power supply adapter from the power supply network.
- Do not exceed the maximum filling volume of the internal chamber indicated in these instructions for use.

## DESCRIPTION OF THE DEVICE

The device is a portable micronized nasal irrigation device for aerosol therapy. The device is composed of a base unit containing a micro pneumatic compressor and a

rechargeable battery. On the upper part, there is the internal chamber with the related nebulization bell.

## GROUPS OF PATIENTS

The device guarantees a highly effective inhalation therapy for patients of all ages, from young children to adults. The device is suitable for home use.

## PACKAGE CONTENTS

The package contains (Fig. A)

- Base unit (1), fitted with power button (2), battery monitoring led (3), micro USB socket (4), micro USB socket cap (10)
- internal chamber (5) fitted with O-ring seal (11)
- Atomiser (6)
- External chamber (7)
- Micro USB cable (8)
- Power supply adapter (9)
- Instructions for use
- Transport case

## PREPARATION AND USE

After removing the device from the packaging, check that it is intact and there is no visible damage.

Before use, proceed with the operations described in the chapter, cleaning, disinfection and sterilisation.

Before use, check that the battery is charging by pressing the power button (3), if the light button flashes, the device will charge, if the light remains lit, the device is ready.

The charge level of the internal battery of the device is indicated by the 3 blue leds placed at the base. The 3 lit leds indicate that the battery is fully charged. If the device is discharged, the blue leds are off and the light button flashes.

When the device is discharged, lower the cap which covers the micro USB socket and connect it to the power supply adapter.

While charging, the three blue leds flash, when the battery is charged, the leds

remain lit. While recharging, the device does not operate and it cannot be used.

After recharging, remove the cable and close the cap of the micro USB socket.

## INSTRUCTIONS FOR USE

*The following operations will be performed with the power supply cable disconnected and the cap of the micro USB socket is well closed.*

- 1) Position the internal chamber (5) on the base unit (1), press it and rotate it clockwise until it is blocked (Fig. F)
- 2) Pour the liquid (maximum 15 ml) in the internal chamber (Fig. B) checking the correspondence of the liquid with the graded scale
- 3) Insert the atomiser (6) into the tube of the tray: a “click” will indicate that the final position has been reached (Fig. G). The reading of the graded scale should only be made when the atomiser is not inserted in the internal chamber.
- 4) Insert the external chamber (7) by turning it slightly clockwise while closing it to cover the sealing gasket.
- 5) Bring the device close to the nose by laying the top outlet on one nostril (Fig. C). During use, sit upright and in a relaxed position.
- 6) To obtain the nebulisation, press the power button (2), to stop the treatment, release the button
- 7) The nebulized liquid comes out from the hole of the external chamber (7) while in the cavity between the external chamber and the internal chamber, the residual liquid is collected (Fig. D)
- 8) During washing, breathe normally with your nose, keeping the opposite nostril blocked with one finger. If necessary, blow from the nostril during treatment to help expel mucus. Alternate nostrils as needed
- 9) The nebulization is very quick and lasts only a few minutes, when the liquid contained in the device is exhausted, release the activation button of the device

## INDICATIONS

### Physiological washing of the nasal cavities

from 0 to 3 years

to solve the inability of the child to blow his/her nose

with saline solution, thermal water or hypertonic solution (Rinowash solution)

at any age	physiological washing and therapeutic adjuvant	with saline solution, thermal water or hypertonic solution (Rinowash solution)
------------	------------------------------------------------	--------------------------------------------------------------------------------

<b>Therapeutic washing</b>		
At any age	in the course of rhinitis, rhinosinusitis, sinusitis, adenoiditis, secretory otitis media and/or as indicated by the general practitioner	with medicines according to the therapy prescribed by the physician

## **CLEANING, DISINFECTION AND STERILISATION**

*The penetration of liquids into the hand piece results in damage to the electrical components and hazards to the safety of the user.*

Before proceeding with the operations indicated below, remove the accessories from the hand piece, to perform these operations, make sure that the power supply cable is not connected to the USB socket of the device.

If it is necessary to clean, disinfect or sterilise this medical device using methods other than those indicated in this chapter, please check the technical data sheet available from the manufacturer.

To avoid any risk of microbial contamination, the device must be cleaned, disinfected and, if necessary, sterilised after each treatment and before use by carefully following the instructions.

Before being cleaned, disinfected or sterilised, the components should be disassembled according to the indications:

- 1) Check that micro USB socket cap (10) of the base unit is well closed
- 2) Turn the internal chamber (5) counter-clockwise, fitted with atomiser and external chamber, and lift it off the base unit.
- 3) Remove the external chamber (7) with a rotation movement
- 4) Remove the atomiser (6) by pushing the flap (Fig. E)

### **BASE UNIT**

*Check that the micro USB socket cap (10) of the base unit is well closed before performing the following operations*

**Cleaning:** it cannot be cleaned under running water  
Use a soft, damp cloth; a neutral, non-abrasive detergent can be used.

**Disinfection:** Use a sponge or a soft cloth moistened with a cold disinfectant solution (such as, for example, a solution with a sodium hypochlorite concentration not exceeding 2%), following the manufacturer's instructions. Ensure that the solution is compatible with the material of which the base unit is made.

**Sterilisation:** the sterilisation cycle is not provided.

### **INTERNAL CHAMBER, ATOMISER AND EXTERNAL CHAMBER**

**Cleaning:** they can be washed with warm water (approx. 40°C) using detergent (dosed according to the detergent manufacturer's indications) and rinsed under running water. When all parts have been cleaned, they can be dried with a soft cloth.

**Disinfection:** after being washed, they can be immersed into a cold disinfectant solution (such as for example a solution with a sodium hypochlorite concentration not exceeding 2%), following the manufacturer's instructions. It is possible to disinfect the accessories by boiling them in water for up to 10 minutes: the parts must not be in direct contact with the bottom of the pot.

**Sterilisation:** they are compatible with the sterilisation cycle in an autoclave at 121°C for at least 15 min. (max. 30 min.) for a maximum of 20 cycles.

## **MAINTENANCE**

The cleaning, disinfection and sterilisation operations of the nasal shower components allow the continuous control of their appearance and can therefore highlight the need for their replacement.

Repeated cleaning, disinfection and sterilisation cycles may alter the colour or aesthetics of the components; this does not affect the operation of the device.

Sterilisation cycles can compromise the characteristics of the nasal shower material; it is therefore advisable to check the appearance and integrity of the components after each sterilisation.

The device must be replaced if damaged.

In case of prolonged non-use, it is recommended to recharge the battery every 6 months to ensure maximum battery efficiency.

## TROUBLESHOOTING

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
No nebulized liquid comes out from the device	Absence of liquid in the internal chamber	Pour the liquid in the internal chamber (5)
	Button not pressed	Press the button (2)
	Battery empty	Recharge the battery
	The device is connected to the power supply network for the recharge	Disconnect the power supply adapter from the device
	The internal chamber, fitted with atomiser and external chamber, is not well fixed to the base unit.	Position the internal chamber (5) on the base unit (1), press it and rotate it clockwise until it is blocked.
	The atomiser is incorrectly inserted into the internal chamber or it is absent.	Check the presence of the atomiser (6) in the internal chamber (5). Remove it and reposition it.
Loss of liquid at the base of the bell	The external chamber (7) is not well inserted and the gasket is not sealing	Insert the external chamber (7) correctly by slightly rotating it during the insertion, checking that the gasket is completely covered
The power button flashes when pressed	The battery is empty	Recharge the battery
	The device is connected to the power supply network for the recharge	Disconnect the power supply adapter from the device


If, after checking the above, the device does not work or in case of other problems or changes in performance, please check the device at an authorised centre.

Do not tamper with or open the base unit of the device. For repair operations, contact only a technical service centre authorised by the manufacturer. Failure to observe the above provisions may compromise the safety of the device.

## TECHNICAL DATA

Power supply	Internal battery with Lithium Polymer 3.7V 1650mAh
Power supply adapter	input 100-240V 50/60Hz 0.5A output 5V DC 2A
Minimum filling volume	3 ml
Maximum filling volume	15 ml
Max. compressor pressure	1 bar
Max. compressor flow	5 lpm
Weight	260 g
Dimensions	H210 x 60 x100 mm
Noise level	58 dbA
Use	30 min ON / 30 min OFF
Battery duration	45 minutes
Service life of the device	400 hours
Aerosol output(*)	13 ml
Aerosol output rate(*)	10 ml/min
MMAD(*)	49 micron

(\*): values measured according to UNI EN 13544 -1. The indicated values refer to the use of physiological solution (0.9% NaCl): they may vary depending on the medicine used. The indicated values do not apply to medicinal products dispensed in suspension or with high viscosity. In that case, the information must be requested from the supplier of the medicine.

Measured with Malvern Spraytec

## DISPOSAL OF THE DEVICE AND ITS ACCESSORIES

The device and/or its accessories, if decommissioned, should be disposed of in compliance with current regulations and environmental protection laws. Where there are no legal requirements, we recommend recycling. The table lists the various types of components of the device.

All materials used are phthalate free and the presence of natural rubber has not been detected.

PART	MATERIAL
Base unit (1)	Waste from electrical and electronic equipment (WEEE)
internal chamber (5)	PC polycarbonate
Atomiser (6)	PC polycarbonate
External chamber (7)	PC polycarbonate + silicone
O-ring seal (11)	Silicone
Micro USB cable (8)	Waste from electrical and electronic equipment (WEEE)
Power supply adapter (9)	Waste from electrical and electronic equipment (WEEE)
Instructions for use	Paper
Transport case	Nylon
Box	Cardboard

## ELECTROMAGNETIC COMPATIBILITY

This device complies with the requirements of CEI EN 60601-1-2: 2015 for the electromagnetic compatibility of medical devices.

Compliance with EMC standards does not guarantee the total immunity of the product; some devices (cell phones, tablets, PC, etc.), when used near medical equipment, can interfere with the operation.

In case of interference with other devices, remove the device.

## OPERATING CONDITIONS

- temperature: from +5°C to +40°C
- relative humidity: from 15% to 93%
- atmospheric pressure: from 700 hPa to 1060 hPa

## STORAGE CONDITIONS

Temperature: -20°C to +60°C

Relative humidity: <85%

Atmospheric pressure: 700 hPa – 1060 hPa

## **SPARE PARTS**

Use only original spare parts available at authorised dealers.

## **UPDATES TO TECHNICAL CHARACTERISTICS**

In order to continuously improve performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are periodically reviewed and improved. Operating instruction manuals are therefore amended to ensure constant compliance with the characteristics of the devices launched to market. If the instruction manual accompanying this device is lost, a copy corresponding to the supplied device may be obtained from the manufacturer by specifying the data shown on the technical data label.

## **COPYRIGHT**

All information contained herein may not be used for purposes other than originally intended.

This manual is owned by Air Liquide Medical Systems S.r.l. and it may not be reproduced, in whole or in part, without the prior written permission of Air Liquide Medical Systems S.r.l. All rights reserved.

## **WARRANTY TERMS AND CONDITIONS**

The manufacturer is liable for the safety, the reliability and the performance of the device if the latter is used according to the instructions, its intended use and if all repairs are carried out at the manufacturer or at a laboratory authorized by the latter.

If the user fails to comply with the Instructions for use and Warnings contained in this manual, the manufacturer's liability and warranty conditions shall be considered void.

The warranty does not cover the accessories and the parts subject to normal wear, such as the battery.

For further information concerning the warranty, please refer to the enclosed warranty coupon.

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