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alvea

 Please read the instructions for use and warnings carefully before using the device (section 1.3).

1 - IMPORTANT INDICATIONS

1.1 - INTENDED USE

Alvea (nebulization ampoule with mouthpiece or mask), can be used to treat bronchial asthma, chronic bronchitis, pulmonary emphysema, cystic fibrosis and other respiratory pathologies by administering medicines in the form of aerosols. Alvea is designed to spray the medicine into particles of a specific size, so that they can be conveyed by the air flow of the compressor to the lower respiratory tract, in the concentration and amount decided by the doctor.

 Medicine administration should be subject to a prescription defining the type of medication, the doses to be administered and the duration of the treatment.

The Alvea ampoule has been designed and manufactured in Italy according to the provisions of Directive 93/42 EEC on Medical Devices (and subsequent updates) and bears the CE mark.

1.2 - DESCRIPTION OF SYMBOLS USED

Caution read the operating instructions  Manufacturer  Lot number 

Keep dry  Free from phthalates  Free from natural latex 

Temperature limits  Pressure limits  Relative humidity limits 

CE marking  Production date  Do not release to the environment 

Recyclable paper packaging  Medical device  Product reference 

1.3 - WARNINGS

- Carefully read the instructions for use and store them in a safe place.
- The device should only be used for low-respiratory tract therapy. Any other use is strictly prohibited.
- Air Liquide Medicaï Systems S.r.l. is not responsible for the equipment that is not manufactured by the same and with which the Alvea ampoule is used. Use devices manufactured in accordance with the provisions of Directive 93/42 EEC on Medical Devices (and subsequent updates) and bearing the CE marking.
- Some parts, due to their small size, could be swallowed causing suffocation.
- It is advisable to check in the table (Chapter 7) if there are any materials due to which some form of allergic reaction has occurred in the past.
-  - Carefully read the operating instructions of the device with which the Alvea ampoule is used to verify its compatibility before using it (see 6 - TECHNICAL DATA).
- Please observe the warnings of the manufacturer of the device used with the Alvea ampoule.
- The device can be used by a child only under the supervision of an adult who is aware of these Instructions.
- Medication in the form of aerosol should be administered only upon medical advice, i.e. according to doctor's prescription defining the type of medicine, the doses to be administered and the duration of the treatment.
- Do not use Alvea while bathing or showering.

1.4 - PATIENT GROUPS

Alvea guarantees highly effective inhalation therapy for patients of any age, from children of paediatric age to adults.

1.5 - PACKAGE CONTENTS

Alvea ampoule is available in four versions:

- Alvea: with adult buccal mask
- Alvea Junior: with paediatric mask
- Alvea Basic: with mouthpiece
- Alvea Family: with paediatric mask, adult buccal mask and mouth piece.

The package contains (Fig. A) the parts listed below.

A nebulizer ampoule consisting of:

- upper body (1) with atomizer (2)
- lower body (3)
- lid (4) with inhalation valve (5) and exhalation valve (6)

A patient interface device, which, depending on the model, can be:

- mouthpiece (7)
- adult mask (8)
- paediatric mask (9)

The lower body (3) is equipped with a terminal to which the air supply hose is connected.

2 - PREPARATION AND USE

Alvea requires connection to an air source which can be the compressed air produced by an aerosol feeding device, medicinal air derived from a centralized system. In all cases it is necessary to ensure consistency of the performance of the air source with the operating data of the ampoule (see section 6 - TECHNICAL DATA).

Read chapter 3 before use.

To use Alvea you must perform the following actions in the order indicated:

- verify that the atomizer (2) is correctly inserted into the upper body of the ampoule (1)
- the lower body of the ampoule (3) is equipped with a graduated scale indicating the amount of liquid inserted: this has only the purpose of giving an approximate indication of the content. For precise dosage use a graduated syringe. Pour the liquid into the tray (Fig. B) according to the amount recommended by your doctor
- hook the upper body of the ampoule (1) onto the lower body (3) by rotating it until you hear a click. (Fig.C-1)
- Check the correct positioning of the valves (5) and (6) in the lid (4)
- hook the lid (4) to the upper body of the ampoule (1) (Fig.C-2)
- Attach the mask (8) (9) or the mouthpiece (7) to the upper body of the ampoule (1) (Fig. C-3)
- Connect the lower terminal of the ampoule to the tube and connect the other end of the tube to the air intake of the aerosol
- put the device into operation
- during treatment the ampoule should be held in a vertical position in order to ensure its proper functioning
- If the operations described above have been carried out correctly, the nebulized liquid will come out from the mask/mouthpiece.

The volumetric buccal mask (8) should be placed gently against the profile of the mouth, while the nose should rest on the outside of the mask. The mouthpiece, on the other hand, (7) should be inserted between the lips: for the effectiveness of the therapy it is necessary that the user inhales through the mouth and exhales through the nose. If the user exhales through the mouth, the valve (6) placed on the lid allows for the exhalation of the exhaled air (Fig. (D)).

The paediatric mask (9) should be placed gently on the patient's face to cover both the nose and the mouth.

Do not clog the valve during use (Fig. (D)).

If the medicine contained in the ampoule has finished or you wish to stop temporarily the treatment, turn off the aerosol supply device or close the air supply. At the end of therapy, a small amount of fluid may remain inside the ampoule; this does not affect the treatment.

3 - CLEANING, DISINFECTION AND STERILIZATION

 To carry out these operations, make sure that the power cord is not connected to the electrical system and that the Alvea ampoule and the hose are not connected to the device.

 If it is necessary to clean, disinfect or sterilize this medical device using methods different than those described in the chapter, check whether this is possible, by consulting the technical data sheet that can be requested to the manufacturer.

To avoid the risk of microbial contamination, Alvea should be cleaned, disinfected and, if necessary, sterilized after each treatment and before use, following the instructions carefully. Before being cleaned, disinfected or sterilised, the parts of the device must be disassembled as follows:

- unplug the hose
- remove the mask or the mouthpiece (7, 8 or 9) from the ampoule (1, 3, 4)
- remove the lid (4) from the upper body of the ampoule (1)
- remove the upper body of the ampoule (1) from the lower body (3)

3.1 - CLEANING

Once the ampoule parts have been disassembled, masks and mouthpiece can be washed with hot water (about 40 °C) using detergent (dosed according to the

manufacturer of the detergent itself) and rinsed under running water. When all the the parts that have been cleaned can be dried with a soft cloth and then stored.

3.2 - DISINFECTION

After washing, all components of the ampoule, masks and mouthpiece can be soaked in a cold disinfectant solution (such as a solution with a sodium hypochlorite concentration above 2%), following the instructions indicated by the manufacturer of the same. It is possible to disinfect Alvea components by boiling them in water for up to 10 minutes: the components must not come in direct contact with the bottom of the pot. The boiling cycle may alter the colour or appearance of the components; this does not affect the operation of the device.

3.3 - STERILIZATION

The components of the ampoule, masks and mouthpiece are compatible with the autoclave sterilization cycle at 121 °C for at least 15 min. (max 30 min.) for up to 20 cycles.

The sterilization cycle is recommended when the device is used in a hospital environment or when it is used by several patients.

4 - MAINTENANCE

Cleaning, disinfection and/or sterilization operations of Alvea components allow for continuous monitoring of their intact state and can therefore give evidence of the need for replacement thereof. Repeated sterilization cycles may affect the characteristics of the material of the ampoule and masks; it is, therefore, recommended to check after each sterilization the appearance and intact state of the components. Alvea must be replaced when its deterioration is obvious.


5 - TROUBLESHOOTING

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
The nebulized medicine does not come out from the ampoule	- Device switched off	- Turn on the appliance
	- The liquid was not poured into the ampoule	- Pour the liquid into the ampoule
	- The liquid was poured in excessive amounts into the ampoule	- Reduce the amount of medicine poured into the ampoule below the maximum capacity
	- The ampoule is not assembled correctly	- Verify that the atomizer (2) is well inserted and attached to the upper body of the ampoule (1). Check that the air comes out from the small hole on the top of the base of the ampoule and screw the upper body (1) into the base (3) properly.
	- The air outlet hole of the base (3) is clogged	- Clean the hole
	- The air supply hose is not connected properly	- Check the correct fitting of the hose
The performance of the ampoule does not comply with the stated	- The air source does not provide the pressure and flow required by the ampoule	- Check the characteristics of the air source
The medicine does not reach the mouth or impaired breathing	- The mask (8) or (9) or the mouthpiece has not been fitted correctly	- Fit the mask or the mouthpiece correctly
	- The valve (Fig. D) is locked	- Unlock the valve


6 - TECHNICAL DATA

Minimum supply pressure	2 bar (200 kPa)
Maximum supply pressure	3 bar (300 kPa)

Air flow at minimum pressure	7 l/min
Air flow at maximum pressure	15 l/min
Maximum ampoule capacity	9 ml
Medicine delivery at minimum pressure (*)	0,20 ml/min
Medicine delivery at maximum pressure (*)	0,45 ml/min
MMAD at minimum pressure (**)	3,6 μm
MMAD at maximum pressure (**)	3 μm

Compressed air source
2 Bar and 7l/min

Compressed air source
3 bar 15 l/min

() NOTE 1: The values indicated refer to the use of saline solution (0.9% NaCl): may vary depending on the medicinal product used. The values indicated do not apply to medicinal products supplied in suspension or with high viscosity. In this case, information should be requested from the supplier of the medicine.*

() NOTE 2: 50% by volume of atomized particles from the ampoule is smaller than the Mass Median Aerodynamic Diameter (MMAD). This value indicates the ability of the ampoule to produce small particles: clinical studies have shown that the efficacy of the ampoule is measured by the ability to nebulize the medicinal product into particles with a diameter of between 1 and 5 microns to allow maximum deposition on the lungs, thus increasing the effectiveness of the treatment.**

The specified MMAD values were detected with the Malvern Spraytec tool. The diagrams provide in order the percentage of the volume of spray solution transported by particles with a diameter of less than the corresponding value in pm reported in abscissa.

7 - DISPOSAL OF ALVEA AMPOLLE

All parts must be disposed of in accordance with current regulations on plastic containers. The table shows the plastics used. All materials used do not contain phthalates and natural latex was not detected.

PART	MATERIAL
Upper body (1)	Polypropylene
Atomizer (2)	Polypropylene
Lower body (3)	Polypropylene
Lid (4)	Polypropylene
Inhalation valve (5)	Silicone
Exhalation valve (6)	Silicone
Mouthpiece (7)	Polypropylene
Adult mask (8)	Silicone
Paediatric mask (9)	Silicone

8 - OPERATING CONDITIONS

- temperature from +5°C to +40°C
- relative humidity from 15% to 93%
- atmospheric pressure from 700 hPa to 1060 hPa

9 - STORAGE CONDITIONS

- temperature from -25 °C to +50 °C
- relative humidity less than 93%
- atmospheric pressure from 700 hPa to 1060 hPa

10 - WARRANTY TERMS AND CONDITIONS

The Alvea ampoule is covered by a 2-year warranty. The manufacturer shall be responsible for the safety, reliability and performance of the device if it is used in accordance with the instructions for use and its intended use and if repairs are carried out by the manufacturer itself or in a laboratory authorised by the same. If the user fails to comply with the operating instructions and the warnings contained in this manual, all liability of the manufacturer and warranty conditions shall be considered void.

11 - UPDATE OF TECHNICAL FEATURES

In order to continuously improve performance, safety and reliability, all medical devices produced by Air Liquide Medical Systems S.r.l. are periodically reviewed and improved.

instruction manuals are therefore modified to ensure their consistency with the characteristics of the devices placed on the market. If the instruction manual accompanying this device is lost, a copy of the version corresponding to the device supplied may be obtained from the manufacturer by specifying the reference data on the technical data sheet.

12 - COPYRIGHT

All information contained in this manual should be used exclusively for the purposes specifically stated. This manual is owned by Air Liquide Medical Systems S.r.l. and it may not be reproduced, either wholly or in part, without the prior written permission of the company. All rights reserved.

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