
4


3

Vitapep Dispositivo per la ginnastica respiratoria
Device for respiratory gymnastic

Vitapep

Dispositivo per la ginnastica respiratoria
Device for respiratory gymnastic

Air Liquide
HEALTHCARE

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ITALIANO

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Symbols

	Warning: read carefully the instructions contained in the manual
	Latex free
	Phthalates free
	Warnings: read carefully the instructions for use
	Temperature limits from -25°C to 50°C
	Relative humidity limits from 30% to 95%
	Atmospheric pressure limits from 760 hPa to 1060 hPa
	Reference code
	Lot number
	Medical device

Warnings

Read carefully the instructions for use and keep them in a safe place.

- There are no contraindications to the use of the device
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established
- The use of the device by a child must always take place under the supervision of an adult who is aware of these instructions.
- The device is personal, for one patient only. The device can be used on the same patient for an indefinite number of treatments. Check after every treatment if the device is complete and not damaged, degradation replace it with a new one.
- The manufacturer has defined a useful life of at least 5 years from the first use.

Functioning principle

Among the single alveoli, the bronchiole and the distal portions of the lung there are collateral canals. Such canals are open under normal conditions, particularly during deep inhaling phases such as sighing which

Data

Vitapep is a medical device for prevention and therapy of secretory problems in patients with acute and chronic respiratory illnesses. Its functioning is based upon a scientifically proven principle: its efficacy is documented in many publications which record the results of clinical tests.

The Vitapep medical device comprises a casing provided with a rotating disk with gauged holes of differing dimensions which can be selected through a suitable key, and a connector provided with a unidirectional valve; the air expired by the patient passes through the selected hole, producing a positive pressure. This positive pressure is obtained only during the expiration phase owing to the connector provided with a unidirectional valve which allows the use of masks of varying types, particularly according to the age of the patient.

Figure 1 shows the various device components.

Figure 4 shows the correct position of the components for the assembly.

Checks to be carried out before using

Before using Vitapep, please ensure that:

- the device is supplied with all its components;
- the device does not show evident signs of damaging;
- the device is clean;
- the selected gauged hole on the rotating disk corresponds to the medical prescription: each one of the 9 selected positions on the rotating disk must be identified by means of a letter (from A to I, in diminishing diameter order); by inserting the special key into the disk radial incisions it is possible to select the prescribed hole by rotating the disk in anti-clockwise manner, as indicated by the arrows. See figure 2. The disk is provided with stopping notches to allow correct positioning of each gauged hole.
- the unidirectional valve must be mounted with the blue side on the body as shown in the figure;
- the mask must adhere perfectly to the patient's face;

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Use of the device

The device, which must be assembled as described above, must be applied to the patient and the mask must adhere perfectly to the face in order to avoid losses of exhaled air between the mask and the face itself. If the mask dimensions are not suitable for the face of the patient, the device (or at least the mask) must be replaced with a suitable one. The treatment consists of deep inspirations and subsequent expirations: after 10-15 breathing cycles, the device must be disconnected from the patient's face who must alternately cough and breath out in a forced manner to expel the secretions and thus free the lungs.

The gauged hole in the specific patient produces during the expiratory phase the necessarily positive pressure for treatment efficacy.

This is carried out by the doctor using a specific Vitapep equipment with a pressure measuring instrument.

For information purposes, the necessary pressure for treatment efficacy is as follows:

- in common cases: 10-15 hPa
- in emphysema cases: 5-8 hPa
- in atelectasis cases: 12-20 hPa

Note: 1 hPa is equivalent to 1 cm H₂O

Use of the device

For treatment frequency, reference to the medical prescription must be made. For information purposes, the recommended frequency is as follows:

- in chronic cases: 10-15 minutes ¼ times per day;
- in acute illness cases: more often, according to secretion quantities;
- in atelectasia and prophylaxis: every hour.

Duration of the treatment:

Treatment duration must also be prescribed by the doctor. For chronic patients, the treatment may last for the entire life, but during the periods in which the secretions problems are not observed, the treatment may be stopped. During post-operative prophylaxis, the treatment must be continued through the entire period in which the patient is immobile.

Patient posture:

The best position from an efficacy point of view of the treatment is when the patient is seated on a chair with his elbows on the table and keeping his torso erect. For bed-ridden patients, use a slanting back stand with arm rest for the elbows.

Note: the treatment will not produce the

expected efficacy is:

- the mask does not adhere perfectly to the face;
- expiration is too forced as this would reduce lungs' volume;
- frequency and duration do not follow the medical prescription.

Maintenance consists of cleaning the components: any non-corrosive detergent product may be used. Components cleaning has to be performed after every use. Vitapep and its component must not have hot or autoclave sterilization.

We recommend disinfecting and cleaning periodically (at least once a week) the disk with the gauged holes taking great care in positioning the "O" gasket correctly in their proper housing (see figure 4); follow the instructions described above for subsequent selection of the prescribed hole.

Disposal

All components must be disposed in plastic containers.

The table shows the plastic materials used. All the materials used do not contain phthalates and the presence of natural latex is not detected.

COMPONENT	MATERIAL
Body with rotating disk support	ABS
Rotating disk	ABS
Fixing screw	Acciaio
Connector with unidirectional valves	ABS and silicone
Fixing OR gasket	Silicone
Key	ABS and steel

Working and storage conditions

- Temperature from -25°C to +50°C
- Relative humidity from 30% to 95%
- Atmospheric pressure from 760 hPa to 1060 hPa

All rights reserved. All information contained in this publication may not be used for purposes other than the original ones. Partial or entire reproduction of this publication is forbidden without the written consent of Air Liquide Medical Systems S.r.l.

With a view of maintaining and improving the standards of production, functions and increased reliability, the apparatuses of Air Liquide Medical Systems S.r.l. are periodically revised and updated. For this reason, the content of this publication are subject to modification without prior notification. It is possible to obtain the updated version of this document should it become obsolete, by contacting the local distributor or the area agent and specifying the edition in your possession.

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