

Respireo

Primo P

Simple & ergonomic, Respireo™ Primo P is a **nasal pillow mask**, silent and easy to use.

Available in Vented and Non Vented versions, it is especially adapted for patients who wish to enjoy freedom and independence in their activities.

It was also designed for patients who are uncomfortable with pressure on the bridge of the nose and/or on the forehead, who suffer from claustrophobia or have a moustache.



Simple & ergonomic

Respireo

Primo P

Adaptable tubing

- 360° mobile ball joint system
- 4 possible attachments (on top of the head, on the sides, or freely downward)



Comfortable and stable headgear

- The angle between the two adjustments provides stability on the head and prevents the upper strap from slipping behind the head
- Two-level adjustment (top and bottom of the head) for a customized fit



Ergonomic cushion

- Wide-edge end caps for an optimal seal
- Easy cushion assembly thanks to the mistake-proofing
- Adapted to everyone's morphology (three sizes available)



Respireo™ Primo P is also available in **Non Vented** version

Silent and diffuse leak



Product reference	Description	Size	Distributor code
472105	Respireo™ Primo P	-	KM219900
472106	Respireo™ Primo P Non Vented	-	KM220000

Spare parts:

472210	Respireo™ Primo P cushion x10	S	KM220700
472211	Respireo™ Primo P cushion x10	M	KM220800
472212	Respireo™ Primo P cushion x10	L	KM220900
472213	Main body assembly Respireo™ Primo P x10	-	KM221000
472214	Tube assembly Respireo™ Primo P Vented x10	-	KM221100
472215	Tube assembly Respireo™ Primo P Non Vented x10	-	KM221200

Manufactured by

AIR LIQUIDE MEDICAL SYSTEMS S.r.l
Via Bisceglie, 66 - 20152 Milano (MI) - ITALY
PLANT: Via dei Prati, 62 - 25073 Bovezzo (BS) - ITALY
Tel. +39 030.201.59.11
Fax +39 030.209.83.29

Distributor Contact

AIR LIQUIDE MEDICAL SYSTEMS S.A.
Parc de Haute Technologie
6, rue Georges Besse
92182 Antony Cedex
France

medicaldevice.airliquide.com

Hotline ALMS: +33 179 51 70 01

Air Liquide
HEALTHCARE

Respireo™ Primo P is a reusable nasal pillow mask with calibrated exhalation orifices (Vented), intended to be used at home or in a hospital setting by a single adult patient (weighing more than 30 kg) who has been prescribed for a non-invasive positive pressure ventilation therapy (NPPV), for example a continuous positive airway pressure therapy (CPAP) or Bilevel positive airway pressure therapy (BiPAP).

Respireo™ Primo P Non Vented is a reusable nasal pillow mask without calibrated exhalation orifices intended to be used at home or in a hospital setting by a single adult patient (weighing more than 30 kg) in controlled non-invasive positive pressure ventilation (NIPPV).

Medical devices class IIa CE0051. Please read the user manual.